

# Beconfident Teeth Whitening LED Booster Light

160298

With Beconfident LED Booster Light, you both improve and prolong the results of your teeth whitening. It is helpful with all Beconfident's teeth whitening kits and has been specially developed for use with the Teeth Whitening Pro Kit as the gel contains Ferrous Gluconate, which is activated by the LED light and makes your teeth whitening even more effective.

The Beconfident LED teeth whitening light has an automatic timer of 15 minutes and can be used with the included mouth tray or together with the mouldable mouth trays that are included in our teeth whitening kits and which are also available as replacement custom fit teeth whitening mouth trays.

Battery is included. Before first use - open battery compartment and remove small plastic disc, which is used to separate the batteries during shipping. Replace batteries and replace battery compartment lid and the LED Booster light is ready to use.



## Technical

---

Country of Origin                      Sweden

## Dimensions

---

Weight                                      0.05 kg

## Package Contents

---

- 1 x Beconfident LED Booster Light
- 1 x Double-sided Mouth Tray