

Beconfident Teeth Whitening Pro Kit

125052

Beconfident Pro Kit with Triple Boost Max Teeth Whitening Gel / X4 uses 3 technologies to whiten teeth up to 7 shades whiter in 10 days. The gel contains the latest PAP formula (ϵ -phthalimido peroxy hexanoic acid) + Sodium Bicarbonate + coco sulphate, which work together for the most effective peroxide free teeth whitening treatment available. An even better result is achieved by also using the Beconfident LED Booster Light to improve the whitening effect by at least one shade - it activates the Ferrous Gluconate ingredient in the gel and improves the whitening results. The teeth whitening kit contains 2 x Beconfident original mouldable mouth trays, a gel applicator and post-whitening treatment. The Protect and Repair Gel (10 ml tube) reduces teeth sensitivity after each treatment and helps to remineralise, strengthen and protect the enamel long term. The Protect and Repair gel contains Hydroxyapatite (liquid calcium) and Sodium monofluorophosphate (920 ppm F⁻). Instructions for use are included.



Technical

Content	10 ml ; 10 ml
Carbamide Peroxide	0% (Zero)
Hydrogen Peroxide	0% (Zero)
CPNP Number	3316120
Country of Origin	Sweden

Ingredients

Glycerin, Sodium Bicarbonate, Phthalimidoperoxycaproic acid, Mentha Piperita Oil, Sodium Coco-Sulfate, Xylitol, Chondrus Crispus Powder, Sorbitol, Aqua, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Flower Oil, Ferrous Gluconate, Limonene, Linalool. REPAIR & PROTECT: Aqua, Glycerin, Sorbitol, Calcium Hydroxyapatite, Carbomer, Tetrasodium Pyrophosphate, Aloe Barbadensis Leaf Juice, Cocos Nucifera Oil, Sodium Coco-Sulfate, Petrolatum, Sodium Monofluorophosphate, Sodium Saccharin, Xanthan Gum, Chamomilla Recutita Flower Extract, Mentha Piperita Oil, Tocopheryl Acetate, Polysorbate 20, Phenoxyethanol, Triethanolamine.

Package Contents

- 1 x Beconfident Teeth Whitening Gel Syringe 10 ml – enough for 10 treatments
- 1 x Beconfident Protect & Repair Teeth Gel 10 ml – enough for 10 treatments
- 2 x Customisable Mouth Trays
- 1 x Gel applicator
- 1 x Instructions for Use

Dimensions

Weight	0.05 kg
--------	---------