

Beconfident Teeth Whitening Dual Foam with Fluoride

131097

Beconfident Dual Foam is a teeth whitening, foaming toothpaste with fluoride. During the treatment, Dual Foam replaces your regular toothpaste. The treatment consists of two steps that both clean and whiten your teeth, a perfect alternative for those who cannot or do not wish to whiten their teeth with the help of teeth whitening gel and mouth trays.

Just brush with Dual Foam for 3 minutes morning and evening. Step 1 Cleaning prepares and cleans the tooth surface for 1 minute followed by Step 2 Whitening for 1 minute. Do not rinse your mouth after brushing, but let the active ingredients work.

Dual Foam is enough for two brushes per day for about three weeks and during this time you should preferably not use regular toothpaste. Dual Foam contains Fluorine, which strengthens the enamel and can be used with both a manual brush and an electric toothbrush. Does not contain abrasives.



Technical

Content 2x50 ml
Carbamide Peroxide 0% (Zero)
Hydrogen Peroxide 0% (Zero)
CPNP Number 2052583
Country of Origin Sweden

Ingredients

Cleaning Foam: Aqua, Alcohol, Glycerin, Polysorbate 20, Poloxamer188, Sodium Lauryl Sulfate, Sodium Saccharin, Sucralose, Poloxamer 407, Mentha Piperita Oil, L-Menthol, Sodium Citrate, Potassium Acesulfame, Sodium Benzoate, Caffeine, Limonene, Sodium Fluoride. Whitening Foam: Aqua, Alcohol, Menta Piperita Oil, Sodium Lauryl Sulfate, Sodium Saccharin, Sucralose, Poloxamer 407, Sodium Coco-sulfate, Disodium EDTA, Sodium Citrate, Potassium Hydroxide, Limonene, Linalool.

Package Contents

- 1 x Cleaning Foam (50 ml)
- 1 x Whitening Foam (50 ml)
- Shade Guide
- Instructions for Use

Dimensions

Weight 0.16 kg

beconfident.com.au 17/11/2025