

# Beconfident Teeth Whitening Pro Kit + LED Booster Light

BCBDL001

Beconfident Pro Kit with Triple Boost Max Teeth whitening Gel / X4 uses 3 technologies to whiten teeth up to 7 shades whiter in 10 days. The gel contains the latest PAP formula ( $\epsilon$ -phthalimido peroxy hexanoic acid) + Sodium Bicarbonate + coco sulphate, which work together for the most effective peroxide free teeth whitening treatment available. The teeth whitening kit contains 2 x Beconfident original mouldable mouth trays, a gel applicator and post-whitening treatment. The Protect and Repair Gel (10 ml tube) reduces teeth sensitivity after each treatment and helps to remineralise, strengthen and protect the enamel long term. The Protect and Repair gel contains Hydroxyapatite (liquid calcium) and Sodium monofluorophosphate (920 ppm F-). Instructions for use are included.

With the Beconfident LED Booster Light, you both improve and prolong the results of your teeth whitening. The teeth whitening LED light has been specially developed for use with the Beconfident Teeth Whitening Pro Kit as the gel contains Ferrous Gluconate, which is activated by the LED light and makes your teeth whitening even more effective.



## Technical

Content	10 ml ; 10 ml
Carbamide Peroxide	0% (Zero)
Hydrogen Peroxide	0% (Zero)
Country of Origin	Sweden

### Ingredients

Glycerin, Sodium Bicarbonate, Phthalimidoperoxyacaproic acid, Mentha Piperita Oil, Sodium Coco-Sulfate, Xylitol, Chondrus Crispus Powder, Sorbitol, Aqua, Aloe Barbadensis Leaf Juice Punica Granatum Extract, Chamomilla Recutita Flower Oil, Ferrous Gluconate, Limonene, Linalool. REPAIR & PROTECT: Aqua, Glycerin, Sorbitol, Calcium Hydroxyapatite, Carbomer, Tetrasodium Pyrophosphate, Aloe Barbadensis Leaf Juice, Cocos Nucifera Oil, Sodium Coco-Sulfate, Petrolatum, Sodium Monofluorophosphate, Sodium Saccharin, Xanthan Gum, Chamomilla Recutita Flower Extract, Mentha Piperita Oil, Tocopheryl Acetate, Polysorbate 20, Phenoxyethanol, Triethanolamine, CI 42090.

## Package Contents

- 1 x Beconfident Teeth Whitening Gel Syringe 10 ml – enough for 10 treatments
- 1 x Beconfident Protect & Repair Teeth Gel 10 ml – enough for 10 treatments
- 2 x Customisable Mouth Trays
- 1 x Gel applicator
- 1 x Instructions for Use
- 1 x Beconfident LED Booster Light
- 1 x Double sided mouth Tray