

## Beconfident Teeth Whitening X1 Start Kit + LED Booster Light + Touch-Up Pen

## BCBDL006

Beconfident X1 Start Kit gives you a peroxide-free teeth whitening for that beautiful white smile without any side effects like pain in the gums or tooth sensitivity. The Start Kit contains everything you need and is perfect even if you haven't whitened your teeth before. It's very easy to use and each treatment only takes 15 minutes. You can easily mould the included mouth trays for an individual fit - just like at the dentist.

Teeth whitening is achieved from the removal of surface stains, including through the oxidation of deep seated organic compounds in the enamel and dentin. Beconfident's method is based on a combination of cleaning agents (Surface Active Agents) and active 32% Sodium Bicarbonate. Unlike hydrogen peroxide, which can cause unpleasant tooth sensitivity and pain in the gums, Beconfident's formula is gentle and free from side effects and whitens teeth up to 4 shades whiter in 10 days.

All that is needed to get your dream smile is included in the package: extra large 10 ml peroxide free whitening gel with fresh melon flavour for 10 treatments, 3pcs (one extra) mouth trays that keeps gel in place, a shade guide and clear instructions.

With the Beconfident LED Booster Light, you both improve and prolong the results of your teeth whitening.

LED Booster Light has an automatic timer of 15 minutes and can be used with the included mouth tray or together with the mouldable mouth tray that is included in our teeth whitening kits.

Battery is included. Before first use - open battery compartment and remove small plastic disc, which is used to separate the batteries during shipping. Replace batteries and replace battery compartment lid and the LED Booster light is ready to use.

The Beconfident Touch-Up Cleaning Pen is suitable for cleaning your teeth on the go or when you need a quick whitening top-up. The teeth whitening pen removes stains from coloured food, tea, coffee, red wine and smoking. The Beconfident pen is applied to the teeth using the inbuilt brush making it the perfect whitening tool to target and remove stains from between teeth. Targeted teeth whitening for a beautifully enhanced whiter smile.



## **Technical**

Content Gel: 10 ml; Pen: 2 ml

Carbamide Peroxide 0% (Zero)
Hydrogen Peroxide 0% (Zero)
Country of Origin Sweden

Ingredients

## Package Contents

- 1 x Teeth Whitening X1 Gel 10 ml enough for 10 treatments
- 3 x customisable and reusable mouth trays
- 1 x Instruction Manual
- 1 x Teeth Whitening Shade Guide (printed on box)
- 1 x Beconfident LED Booster Light

beconfident.com.au 18/8/2025

X1 Start Kit: Glycerin, Sodium Bicarbonate, Aroma, Xylitol, Chondrus Crispus Powder, Polysorbate 20, Sorbitol, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Extract, Sodium flouride, Linalool, Limonene. Touch-Up Pen: Glycerin, Sodium Bicarbonate, Aroma, Xylitol, Chondrus Crispus Powder, Polysorbate 20, Sorbitol, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Extract, Sodium flouride, Linalool, Limonene.

- 1 x Double-sided Mouth Tray
- 1 x Touch-Up Pen 2 ml
- 1 x Instruction Manual
- 1 x Teeth Whitening Shade Guide (printed on box) on box lid

beconfident.com.au 18/8/2025