

## Beconfident Teeth Whitening Pro Kit + LED Booster Light + Touch-Up Pen

BCBDL007

Beconfident Pro Kit with Triple Boost Max Teeth whitening Gel uses 3 technologies to whiten teeth up to 7 shades whiter in 10 days. The gel contains PAP (ε-phthalimido peroxy hexanoic acid) + Sodium Bicarbonate + coco sulphate, which work together for the most effective peroxide free teeth whitening treatment available. The teeth whitening kit contains 2 x Beconfident original mouldable mouth trays, a gel applicator and post-whitening treatment. The Protect and Repair Gel (10 ml tube) reduces teeth sensitivity after each treatment and helps to remineralise, strengthen and protect the enamel long term. The Protect and Repair gel contains Hydroxyaptite (liquid calcium) and Sodium monofluorphosphate (920 ppm F-). Instructions for use are included.

With the Beconfident LED Booster Light, you both improve and prolong the results of your teeth whitening. The teeth whitening LED light has been specially developed for use with the Beconfident Teeth Whitening Pro Kit as the gel contains Ferrous Gluconate, which is activated by the LED light and makes your teeth whitening even more effective.

With the Beconfident LED Booster Light, you both improve and prolong the results of your teeth whitening. The teeth whitening LED light has been specially developed for use with the Beconfident Teeth Whitening Pro Kit as the gel contains Ferrous Gluconate, which is activated by the LED light and makes your teeth whitening even more effective.

The Beconfident Touch-Up Cleaning Pen is suitable for cleaning your teeth on the go or when you need a quick whitening top-up. The whitening pen removes stains from many culprits including coloured food, drink and smoking. The Beconfident pen is applied to the teeth using the brush making it the perfect whitening tool to target and remove stains from between teeth. Beconfident peroxide free teeth whitening pen is perfect for removing stains from between teeth, helping to maintain whitened teeth after using the Beconfident whitening kit or to use alongside Beconfident whitening treatments. Targeted teeth whitening for a beautifully enhanced whiter smile.



## **Technical**

Content Gel: 10 ml; 10 ml; Pen: 2 ml

Carbamide Peroxide 0% (Zero)
Hydrogen Peroxide 0% (Zero)
Country of Origin Sweden

Ingredients

Pro Kit: Glycerin, Sodium Bicarbonate, Phthalimidoperoxycaproic acid, Mentha Piperita Oil, Sodium Coco-Sulfate, Xylitol, Chondrus Crispus Powder, Sorbitol, Aqua, Aloe Barbadensis Leaf Juice Punica Granatum Extract, Chamomilla Recutita Flower Oil, Ferrous Gluconate,

## Package Contents

- 1 x Beconfident Teeth Whitening Gel Syringe 10 ml enough for 10 treatments
- 1 x Beconfident Protect & Repair Teeth Gel 10 ml enough for 10 treatments
- 2 x Customisable Mouth Trays
- 1 x Gel applicator
- 1 x Instructions for Use
- 1 x Beconfident LED Booster Light
- 1 x Double-sided Mouth Tray
- 1 x Touch-Up Pen 2 ml

beconfident.com.au 12/7/2025

Limonene, Linalool. REPAIR & PROTECT: Aqua, Glycerin, Sorbitol, Calcium Hydroxyapatite, Carbomer, Tetrasodium Pyrophosphate, Aloe Barbadensis Leaf Juice, Cocos Nucifera Oil, Sodium Coco-Sulfate, Petrolatum, Sodium Monofluorophosphate, Sodium Saccharin, Xanthan Gum, Chamomilla Recutita Flower Extract, Mentha Piperita Oil, Tocopheryl Acetate, Polysorbate 20, Phenoxyethanol, Triethanolamine. Touch-Up Pen: Glycerin, Sodium Bicarbonate, Aroma, Xylitol, Chondrus Crispus Powder, Polysorbate 20, Sorbitol, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Extract, Sodium flouride, Linalool, Limonene.

- 1 x Instruction Manual
- 1 x Teeth Whitening Shade Guide (printed on box) on box lid

beconfident.com.au 12/7/2025