

Beconfident Value Teeth Whitening Kit

BCBDL013

Beconfident's patented Dual-Boost X3 teeth whitening gel uses three different modes of action to brighten teeth up to 6 shades whiter in 14 days. In the first active phase, superficial discolourations are gently removed from the teeth. In the second active phase, a highly efficient whitening ingredient complex removes deeper discolouration from the teeth and visibly brightens them. Following the gel treatment, the additional mint flavoured Tooth Gloss protects teeth from new discolouration occurring and gives them fantastic shine. The result - bright white teeth for a beautiful attractive smile. The package contains 1 x 8 ml teeth whitening gel syringe which is enough for 8 days use, plus 1 x Tooth Gloss Mint (10 ml tube). Instructions for use also included.

Beconfident custom fit Teeth Whitening Mouth Trays can be used together with all Beconfident's different teeth whitening gels. A suitable mouth tray is a prerequisite for a successful teeth whitening result as it is critical that the gel has as optimal contact as possible with the tooth surface. The mouth trays come with detailed instructions on how to easily shape each one to your own upper and lower tooth row. The package contains 2 x Teeth Whitening mouth trays plus gel applicator for treatment of individual tooth surfaces and instructions for use.





Technical

Content Gel: 8 ml; Gloss Mint: 10 ml

Carbamide Peroxide 0% (Zero)

Hydrogen Peroxide 0.1%

CPNP Number 2803188

Country of Origin USA

Ingredients

SYRINGE: Glycerin, Sodium Bicarbonate, Mentha Piperita Oil, Sodium Coco-Sulfate, Xylitol, Chondrus Crispus Powder, Sorbitol, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Sodium Flouride, Chamomilla Recutita Extract, Limonene, Linalool, Urea Peroxide, CI 42090 TOOTHGLOSS: Aqua, Poloxamer 407, Alcohol, Carbopol, Sucralose, Sodium Saccharin, L-Menthol, Sodium Benzoate, Aroma

Package Contents

- 1 x Teeth Whitening Gel Syringe 8ml
- 1 x Tooth Gloss mint 10ml
- 1 x Instruction Manual
- 2 x Customisable Mouth Trays
- 1 x Gel Applicator
- 1 x Instruction Manual

Dimensions

Weight 0.50 kg

beconfident.com.au 15/12/2025