

SIMPLESMILE® Teeth Whitening Start Kit

SS1200

The teeth whitening gel in the Start Kit is based on a technology with an active whitening ingredient (32% bicarbonate) which together with the cleaning substances removes discolourations and gives whiter teeth.

The included mouth trays are adjusted to each teeth row and provide optimal contact between the teeth whitening gel and the teeth for best effect.

The treatment takes only 15 minutes and the provided gel is enough for about 10 applications. We recommend that you complete the treatment until the gel is finished and then repeat when needed, our recommendation is to do a full treatment approximately every three months. The mouth trays can be cleaned and saved for the next time.

Simplesmile® is a safe and well-proven product that complies with all rules and laws regarding teeth whitening products. Tested and approved by dentists.



Technical

Content	10 ml
Carbamide Peroxide	0% (Zero)
Hydrogen Peroxide	0% (Zero)
Country of Origin	Sweden

Ingredients

Gel: Glycerin, Sodium Bicarbonate, Aroma, Xylitol, Chondrus Crispus Powder, Polysorbate 20, Sorbitol, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Extract, Sodium fluoride, Linalool, Limonene.

Package Contents

- 1 x Teeth Whitening X1 Gel 10 ml – enough for 10 treatments
- 3 x Customisable and reusable mouth trays
- 1 x Instruction Manual
- 1 x Teeth Whitening Shade Guide (printed on box) (on Box)

Dimensions

Weight	0.04 kg
--------	---------