

SIMPLESMILE® Touch-Up Pen

SS1202

Easy to bring with you, gentle on teeth and gums and can be applied as often as you wish – anytime and anywhere. A perfect product for maintaining the whiteness of the teeth between treatments.

You simply brush the bleaching gel evenly over the tooth surface to improve the result. It's also possible to bleach one or more teeth a little extra with the help of the pen.

The teeth whitening gel in the Touch-Up pen is based on the same technology as in the Start Kit with the active ingredient 32% bicarbonate. The gel also contains Aloe Vera which is soothing for the gums, as well as Fluorine to counteract caries.



Technical

Content 2 ml
Carbamide Peroxide 0% (Zero)

Hydrogen Peroxide 0% (Zero)
Country of Origin Sweden

Ingredients

Pen: Glycerin, Sodium Bicarbonate, Aroma, Xylitol, Chondrus Crispus Powder, Polysorbate 20, Sorbitol, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Extract, Sodium flouride, Linalool, Limonene.

Package Contents

- 1 x Touch-Up Pen 2 ml
- 1 x Instructions for Use
- 1 x Teeth Whitening Shade Guide (printed on box) (on Box)

Dimensions

Weight 0.02 kg

beconfident.com.au 18/8/2025