

## SIMPLESMILE® Start Kit + Refill + Touch-Up Pen

SSBDL001

The teeth whitening gel in the Start Kit and the Refill pack is based on a technology with an active whitening ingredient (32% bicarbonate) which together with the cleaning substances removes discolourations and gives whiter teeth.

The included mouth trays are adjusted to each teeth row and provide optimal contact between the teeth whitening gel and the teeth for best effect.

The treatment takes only 15 minutes and the provided gel is enough for about 20 applications (10 in Start Kit + 10 in Refill). We recommend that you complete the treatment until the gel is finished and then repeat when needed, our recommendation is to do a full treatment approximately every three months. The mouth trays can be cleaned and saved for the next time.

The package contains teeth whitening gel with a fresh taste of melon (10 ml), 3 mouldable mouth trays, a shade guide (printed on the box) and instructions for use.

Simplesmile® is a safe and well-proven product that complies with all rules and laws regarding teeth whitening products. Tested and approved by dentists.

The Touch-Up Pen is easy to bring with you, gentle on teeth and gums and can be applied as often as you wish – anytime and anywhere. A perfect product for maintaining the whiteness of the teeth between treatments.

You simply brush the bleaching gel evenly over the tooth surface to improve the result. It's also possible to bleach one or more teeth a little extra with the help of the pen.

The teeth whitening gel in the Touch-Up pen is based on the same technology as in the Start Kit with the active ingredient 32% bicarbonate. The gel also contains Aloe Vera which is soothing for the gums, as well as Fluorine to counteract caries.



## **Technical**

Country of Origin

Sweden

Ingredients

Gel: Glycerin, Sodium Bicarbonate, Aroma, Xylitol, Chondrus Crispus Powder, Polysorbate 20, Sorbitol, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Extract, Sodium flouride, Linalool, Limonene. Pen: Glycerin, Sodium Bicarbonate, Aroma, Xylitol, Chondrus Crispus Powder, Polysorbate 20, Sorbitol, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Extract, Sodium flouride, Linalool, Limonene.

## Package Contents

- 2 x Teeth Whitening X1 Gel 10 ml enough for 20 treatments
- 3 x Customisable and reusable mouth trays
- 2 x Instruction Manual
- 1 x Touch-Up Pen 2 ml
- 1 x Instructions for Use
- 3 x Teeth Whitening Shade Guide (printed on box) (on Box)

beconfident.com.au 18/8/2025