

SIMPLESMILE® Start Kit + Refill

SSBDL007

The teeth whitening gel in the Start Kit and the Refill pack is based on a technology with an active whitening ingredient (32% bicarbonate) which together with the cleaning substances removes discolourations and gives whiter teeth.

The included mouth trays are adjusted to each teeth row and provide optimal contact between the teeth whitening gel and the teeth for best effect.

The treatment takes only 15 minutes and the provided gel is enough for about 20 applications (10 in Start Kit + 10 in Refill). We recommend that you complete the treatment until the gel is finished and then repeat when needed, our recommendation is to do a full treatment approximately every three months. The mouth trays can be cleaned and saved for the next time.

The package contains teeth whitening gel with a fresh taste of melon (10 ml), 3 mouldable mouth trays, a shade guide and instructions for use.

Simplesmile® is a safe and well-proven product that complies with all rules and laws regarding teeth whitening products. Tested and approved by dentists.



Technical

Country of Origin

Sweden

Ingredients

Gel: Glycerin, Sodium Bicarbonate, Aroma, Xylitol, Chondrus Crispus Powder, Polysorbate 20, Sorbitol, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Extract, Sodium flouride, Linalool, Limonene.

Package Contents

- 2 x Teeth Whitening X1 Gel 10 ml – enough for 20 treatments
- 3 x Customisable and reusable mouth trays
- 2 x Instruction Manual
- 2 x Teeth Whitening Shade Guide (printed on box) (on Box)